

How to do a rapid COVID-19 self-test

Your step by step guide



 Before you begin, check you have a sealed swab, a registration card and a barcode




 Use the mirror to look at your tonsils (or where they used to be), then gently blow your nose

Dispose the tissue in the bin provided




 Sanitise your hands then open your swab packet, being careful not touch the fabric tip


Ask a helper for scissors if needed



 Rub the swab tip firmly against your tonsils (or where they used to be) 4 times on each side and remove without touching your cheeks, lips or gums



5) Put swab tip into 1 nostril until you feel slight resistance and roll the swab firmly around the inside of the nostril, making 10 complete circles.



 Place the swab (fabric tip first) directly into the tube on the test tube rack in front of you



 Attach your barcode to the registration card and register using your device or the device provided

Scan or type in barcode and skip NHS number if you do not know it



 Dispose of all rubbish and clean entire booth, including mirror, table, device, bin, seat and scissors (if used)



 Put on your face covering, sanitise your hands and leave test area, taking your registration card with you